# HIGH FLOW BIO-X Air & Gas

Filter Cartridges





HIGH FLOW BIO-X sterile gas filters combine proven depth filter technology and a pleated construction to provide retention down to 0.01 micron in gas.

Flow rates typically 2-3 times that of membrane filters make HIGH FLOW BIO-X the filter that can dramatically reduce cartridge usage and installation size within the fermentation, food and beverage industries.

The specially developed PTFE impregnation process imparts greater strength and permanent hydrophobicity to the borosilicate microfibre media. This leads to excellent performance in applications such as the provision of sterile gas in filling machines.


# **Features**

- I High flowing hydrophic PTFE impregnated media
- I Fully validated by aerosolized bacterial and viral challenge
- I Stainless steel inner core
- 100% integrity testable by Valairdata 3 aerosol challenge

# **Benefits**


- I Reduce system size and reduced total cost of ownership.
- I Provides complete process security
- I Strong and robust for extended service life
- I Guaranteed performance in-situ

# Performance Characteristics



# Filtration Stage

### Sterile Gas and Vent Filtration



 I Outer Protection Cage:
 Polypropylene

 I End Caps:
 Polypropylene

 I End Cap Insert:
 Polyethersulphone

I Standard o-rings/gaskets: Silicone

# Food Contact Compliance

Parker domnick hunter's range of HIGH FLOW BIO-X filters are intended for indirect food contact and as such are manufactured from materials suitable for the sterilization of compressed gasses within food and beverage applications. Materials conform to the relevant requirements of the United States FDA 21 CFR part 177 and USP Plastics Class VI – 121°C.

#### **Recommended Operating Conditions**

The maximum differential pressure in direction of flow (outside to in) is 3.5 barg (50.76 psig) at 70 °C (158 °F).

The maximum recommended continuous operating temperature is 70 °C (158 °F).

#### Effective Filtration Area (EFA)

10" (250 mm) Up to 0.38 m<sup>2</sup> (4.09 ft<sup>2</sup>)

#### Sterilization

HIGH FLOW BIO-X cartridges can be in-situ steam sterilized or autoclaved up to 142 °C (287.6 °F) for a maximum of 150 steam cycles.

For detailed operational procedures and advice on cleaning and sterilization, please contact the Technical Support Group through your usual Parker domnick hunter contact.

#### **Retention Characteristics**

The HIGH FLOW BIO-X range of cartridges has been fully validated by aerosol bacterial challenge with challenge levels of 10<sup>12</sup> *Brevundimonas diminuta* per 10" [250 mm] filter cartridge. Independent test work also shows full retention to *MS-2 Coliphage*.


## Integrity Test Data

All cartridges are integrity tested prior to despatch by the aerosol challenge test method using the Parker domnick hunter VALAIRDATA 3.

#### Manufacturing Traceability

Each filter cartridge displays the product name, product code and lot number. Additionally, each module displays a unique serial number providing full manufacturing traceability.

# Ordering information



HBA & HPG HOUSING RANGE AVAILABLE







# WE ARE YOUR 'ONE STOP SHOP' FOR ALL OF YOUR FILTRATION REQUIREMENTS

## LIQUID FILTRATION

- CARTRIDGES
- FILTER BAGS
- VESSELS
- STRAINERS

## **AIR FILTRATION**

- VENT FILTERS
- STERILE AIR
- VACCUM
- AIR-CONDITIONING
- FILTER MEDIA

## **DOMESTIC**

- ULTRAVIOLET
- TREATMENT OF BOREHOLES, SPRINGS & WELLS
- WATER TESTING

## **FUEL FILTRATION**

• DIESEL, AVIATION



Prosep Filter Systems Ltd Unit G19, River Bank Way, Lowfields Business Park, Elland, West Yorkshire HX5 9DN Phone: 01422 377367 email: sales@prosep.co.uk www.prosep.co.uk





