

Proflow™ II Pharmaceutical LBR

Pleated PTFE Membrane Cartridge for
Pharmaceutical Solvent Sterilization

The Proflow II Pharmaceutical Liquid Bacterial Retentive (LBR) cartridge is designed for those who require sterile pharmaceutical fluids and is validated to provide sterility through the retention of *Brevundimonas diminuta* per HIMA guidelines. Along with its sterilizing-grade performance, this cartridge exhibits high air and fluid flow rates due to the 9.3ft² of expanded PTFE membrane it contains.

Proflow II-LBR cartridges are manufactured in a certified Class 10,000 cleanroom where they are flushed with non-pyrogenic deionized water and 100% integrity tested as well as serialized to provide complete product traceability.



Contact Information

Parker-Hannifin Corporation
domnick hunter
Process Filtration - North America
2340 Eastman Avenue
Oxnard, California, USA 93030

Toll Free: +1 877 784 2234
Phone: +1 805 604 3400
Fax: +1 805 604 3401
dhpsales.na@parker.com

www.parker.com/processfiltration

Benefits

- Fully validated for sterilizing-grade performance
- Broad chemical compatibility allows use with most pharmaceutical solvents
- Long service life minimizes change-out frequency
- Low differential pressure reduces energy costs
- Integrity tested to ensure quality
- Serialized for immediate traceability

Applications

- Solvent sterilization
- Extraction/crystallization solvents
- Hot WFI filtration
- Aseptic packaging air / sterile air supply
- Fermenter/tank vents



ENGINEERING YOUR SUCCESS.

Proflow™ II Pharmaceutical LBR

Pleated PTFE Membrane Cartridge for Pharmaceutical Solvent Sterilization

Liquid Flow Rates*

10 inch (250mm) cartridges
Liquid: 1.9gpm/psid
(10.43 lpm/100 mbar)

5 inch (130mm) cartridges
Liquid: 1.0gpm/psid
(5.5 lpm/100 mbar)

Materials of Construction

Membrane: PTFE
Support layers: Polypropylene
Structure: Polypropylene

All components meet USP XXIII Class VI-121°C criteria and are thermally bonded to assure integrity and minimize extractables.

Effective Filtration Area

10 inch (250mm) cartridges
9.3 ft² (0.86 mm)

5 inch (130mm) cartridges
4.6 ft² (0.43 mm)

Bacteria Retention

Proflow II LBR cartridges will provide a sterile effluent when challenged with greater than 10⁷ CFU of *Brevundimonas diminuta* per cm² of EFA in accordance with the test methodology recommended by HIMA and defined by ASTM method F838-05.


Steam Sterilization

Cartridges may be steamed or autoclaved for over 20 one-hour cycles @ 275°F (145°C). Cartridges may also be chemically sanitized in place using common sanitizing agents.

Maximum Differential Pressure

Forward:
80psid (5.5 bar) @ 75°F (24°C)
40psid (2.8 bar) @ 180°F (82°C)

Reverse:
50psid (3.4 bar) @ 75°F (24°C)



Ordering Information

Each cartridge is identified with a product number, pore size, and lot number for traceability.

34 — [] 0 [] [] [] — 002 — [] — L

Insert Style	
CODE	DESCRIPTION
5	Encapsulated 316L Stainless Steel
6	Encapsulated Polysulfone (Standard)

End Fitting	
CODE	DESCRIPTION
2	226/Flat
3	222/Flat
7	226/Fin
8	222/Fin

Nominal Length	
CODE	FITTING
05	5" (130mm)
10	10" (250mm)
20	20" (500mm)
30	30" (750mm)
40	40" (1000mm)

O-Rings	
CODE	MATERIAL
1	EPDM
2	Silicone (recommended)
4	Viton®
5	Encapsulated Viton
6	Encapsulated Silicone
N	None

Proflow is a registered trademark of Parker-Hannifin Corporation.
Viton is a registered trademark of E.I. duPont de Nemours & Co. Inc.

© 2008 Parker -annifin Corporation
domnick hunter Process Filtration - North America
All Rights Reserved
DS_LS_Proflow II LBR Cart. 02/12 Rev. B



ENGINEERING YOUR SUCCESS.

PROSEP FILTER SYSTEMS LTD
Unit G19
River Bank Way
Lowfields Business Park
Elland
West Yorkshire
HX5 9DN


Tel: 01422 377367

Fax: 01422 377369

Email: enquiries@prosep.co.uk

www.prosep.co.uk

Map and Directions to Prosep Filters Limited



Leave M62 at Junction 24.

At roundabout adjacent to Cedar Court Hotel take 2nd exit onto dual carriageway (A629), signposted Halifax.

Take 1st exit slip road.

At roundabout at end of sliproad, take 3rd exit off.

This is the entrance to Lowfields Business Park.

Proceed straight over 1st roundabout.

At next roundabout take 2nd exit onto River Bank Way - Prosep Filters can be found on the left after the S-bend.

[Link to Google Maps](#)